

PD Drainage Bags



DB-301 (3.5L Drainage Bag)

- 3.5 Litre Drain Bag
- Input port for secure connection of the drain line spike of the Y administration set
- outlet port with clamp to empty the bag
- Calibration marking for measurement of the drain volume
- Hooks and string for hanging the drain bag



DB-302 (Neonatal Drainage Bag)

- 500 ml Pediatric Meter with 2L overflow bag
- Input port for secure connection of the drain line spike of the pediatric Y administration set
- outlet valve to empty the meter
- Outlet tubing to empty the overflow bag
- Calibration marking for measurement of the drain volume
- Hooks and string for hanging the drain bag



DB-1401 - 15L Drainage Bag

- 15 L Drainage Bag with inflow and outflow
- Input spike for secure connection to the cyclor drain line and other similar ports
- Outlet tubing with clamp to empty the drain bag



DC-D302D - Dual 15L Drainage Bag Set

- 2- 15 L Drainage Bags
- Input spike for secure connection to the cyclor drain line and other similar ports

