






Hemodialysis

Double Lumen Catheters



- Made of polyurethane for outstanding biocompatibility
- Soft tip reduces trauma to vessel upon introduction
- Radiopaque catheter and tip facilitates vessel location
- Soft hub allows for increased patient comfort
- Catheter kit includes catheter, dilator, guidewire, introducer needle and injection caps

Straight Extensions

DBLK-12-13 (Size 12Fr x 13cm)
DBLK-12-15 (Size 12Fr x 15cm)
DBLK-12-20 (Size 12Fr x 20cm)

Curved Extensions

DBLXK-12-13 (Size 12Fr x 13cm)
DBLXK-12-15 (Size 12Fr x 15cm)
DBLXK-12-20 (Size 12Fr x 20cm)



Medionics International Inc.

114 Anderson Avenue, Markham, Ontario L6E 1A5 Canada

Tel: (905) 472-6544 Fax: (905) 472-6549

info@medionics.com www.medionics.com



Made in Canada