


# Pediatric IPD Set

## PDS-1100PLS



### **Specifications:**

- ***Y-type pediatric IPD set with 150 ml vented Burette***
- ***For use with dialysis solution bag equipped with luer lock input port***
- ***Drip chamber for viewing drain flow rates***
- ***Rotating Male Luer lock patient line connector***
- ***Roller clamps with super fine flow control***
- ***Drain spike for attaching to drain bag with spike port, such as Medionics' DB-301, DB-302***



Medionics International Inc., 114 Anderson Avenue, Markham, Ont L6E 1A5 Canada  
Tel: (905) 472-6544 Fax: (905) 472-6549 info@medionics.com [www.medionics.com](http://www.medionics.com)



Made in Canada  
Rev 1-03/28/2016