


Pediatric IPD Set

PDS-1200BCCH



Specifications:

- **Y-type pediatric IPD set with 2 – 150 ml vented Burettes**
- **For use with two dialysis solution bags equipped with luer lock input ports**
- **Drip chambers for viewing fill and drain flow rates**
- **Rotating Male Luer lock patient line connector**
- **Roller clamps with super fine flow control**
- **Drain spike for attaching to drain bag with spike port, such as Medionics' DB-301, DB-302**



Medionics International Inc., 114 Anderson Avenue, Markham, Ont L6E 1A5 Canada
Tel: (905) 472-6544 Fax: (905) 472-6549 info@medionics.com www.medionics.com



Made in Canada
Rev 1-03/28/2016