


PD CATHETER KIT



Specifications:

- *Kit versions are available for all catheter types*
- *Kit is used to insert catheters percutaneously using modified Seldinger technique*
- *It contains one peritoneal catheter, one plastic catheter adapter, one catheter clamp, one sealing cap, one Pull-Apart Sheath/Dilator, one 18ga introducer needle, one J/Straight 0.038" Guidewire, one peritoneal tunneling stylet, one #11 scalpel, 4 pcs of gauze sponge and one 10cc syringe*
- *Kit is designed to insert Medionics' 4.95 mm OD catheters*



Medionics International Inc., 114 Anderson Avenue, Markham, Ontario L6E 1A5 Canada
Tel: (905) 472-6544 Fax: (905) 472-6549 info@medionics.com www.medionics.com



Made in Canada