

# ***PD TEMP***

## **Curved Stylet Catheter for Acute Peritoneal Dialysis**



**RCT-288CRK (Size 2.5 x 3.5 x 280mm, Perforations 8cm)**



### ***Specifications:***

- ***Quick, easy and safe***
- ***Made of Polyamide for outstanding tissue compatibility***
- ***Precision hole punching without burrs to prevent catheter occlusion***
- ***Bent Abdomen segment helps to advance catheter in the deep pelvic area***
- ***Elbow connector prevents catheter kinking***
- ***Catheter marking to show direction of distal end to facilitate insertion***
- ***Sampling port for injections during treatment***



Medionics International Inc., 114 Anderson Avenue, Markham, Ontario L6E 1A5 Canada  
Tel: (905) 472-6544 Fax: (905) 472-6549 info@medionics.com [www.medionics.com](http://www.medionics.com)



**Made in Canada**