


Introducer Set

for

Percutaneous Insertion of PD Catheters



Specifications:

- *Introducer Set is used to insert catheters percutaneously using modified Seldinger technique;*
- *It contains one Pull-Apart Sheath/Dilator, one 18ga introducer needle, one J/Straight 0.038" Guidewire and one 10cc syringe;*
- *Designed to insert Medionics' 4.95 mm OD peritoneal catheters;*
- *Available for adult, pediatric and neonatal catheters*



Medionics International Inc., 114 Anderson Avenue, Markham, Ontario L6E 1A5 Canada
Tel: (905) 472-6544 Fax: (905) 472-6549 info@medionics.com www.medionics.com



Made in Canada