

CATHETER ADAPTER

*& accessories for
PD Catheters*

Plastic Catheter Adapter PCA-1010



Injection Cap ICL-100



Catheter Clamp PCA-1030



Specifications:

- ***Plastic Catheter Adapter is an inexpensive alternative to the Titanium Adapter. It allows peritoneal catheters to connect securely to CAPD administration sets***
- ***Catheter Clamp allows easy clamping of the peritoneal catheter***
- ***Injection Cap is used to cap off the catheter adapter. It is self sealing and provides for administering medication***



Medionics International Inc., 114 Anderson Avenue, Markham, Ontario L6E 1A5 Canada
Tel: (905) 472-6544 Fax: (905) 472-6549 info@medionics.com www.medionics.com

