

reliant™

Reinforced PD Catheters

Reliant™ reinforced catheters are Medionics' next generation of PD catheters uniquely designed with additional mass of silicone embedded into the intra-abdominal segment to keep the catheter low in pelvic gutter. Reinforcement helps to prevent catheter migration and stabilizes PD access.

Silicone reinforcement of the preperitoneal cuff strengthens cuff anchorage into the abdominal wall, prevents leakage and helps maintain catheter tip in the true pelvis.

Double layered up/down structure minimizes effects of bowel and omental adhesions. A secondary set of recessed holes provide unimpeded flow even when catheter migration does occur.

Order No. for Reliant™ catheter,
add Suffix "RR" to catheter catalog number.




**Tenckhoff
Reliant™ Catheters**

Order No. for Reliant™ catheter,
add Suffix "RR" to catheter catalog number.



**Tenckhoff Coil
Reliant™ Catheters**

Order No. for Reliant™ catheter,
add Suffix "RR" to catheter catalog number.



**Swan Neck
Reliant™ Catheters**

Order No. for Reliant™ catheter,
add Suffix "RR" to catheter catalog number.



**Swan Neck Coil
Reliant™ Catheters**