


Straightening Stylets

for

PD CATHETERS



Specifications:

- *Stylet is used to straighten PD catheters during catheter insertion*
- *It is available in different lengths for use with catheters of different lengths*
- *Made of stainless steel*
- *Designed to fit Medionics' PD Catheters*
- *Sterile, autoclavable for possible reuse*



Medionics International Inc., 114 Anderson Avenue, Markham, Ontario L6E 1A5 Canada
Tel: (905) 472-6544 Fax: (905) 472-6549 info@medionics.com www.medionics.com

