


Swan Neck PD Catheters

Adult or Pediatric



Specifications:

- Medionics' Swan Neck catheters are designed with a bent segment, adjacent to the cuffs, in its natural, unstressed condition
- The bent segment permits the catheter to be mounted in a tunnel formed through abdominal wall in relatively unstressed configuration, with both ends pointing downwardly resulting in reduced tunnel infection
- Medionics' swan neck catheters are designed to form a substantial angle in its bent segment to avoid the danger of kinking the catheter
- Medionics' swan neck catheters are designed to provide optimum flow for improved dialysis clearances

