

TROCAR

For

PD Catheter implantation



Specifications:

- *Tenckhoff Trocar is specially designed for insertion of cuffed PD catheters*
- *It is a tapered bivalve instrument with a handle on each half and an obturator with a Trocar point*
- *Designed to insert Medionics' 4.95 mm OD peritoneal catheters*
- *Sterile, autoclavable*



Medionics International Inc., 114 Anderson Avenue, Markham, Ontario L6E 1A5 Canada
Tel: (905) 472-6544 Fax: (905) 472-6549 info@medionics.com www.medionics.com

