


# TITANIUM ADAPTER

for  
PD Catheters


Catalog Number TA-200



Assembled TA-200

Lock Nut

Main Body



## Specifications:

- 2 Part Peritoneal Dialysis Catheter Adapter with lock nut
- It is used to terminate the external portion of peritoneal catheters
- It allows peritoneal catheters to connect securely to CAPD administration sets
- Designed to fit Medionics' 4.95 mm OD PD Catheters
- Lock nut secures catheter to the Titanium Adapter main body
- Manufactured with titanium



Medionics International Inc., 114 Anderson Avenue, Markham, Ontario L6E 1A5 Canada  
Tel: (905) 472-6544 Fax: (905) 472-6549 info@medionics.com [www.medionics.com](http://www.medionics.com)

