


# MY SYSTEM

## A Non-wearable CAPD "Y" System

*Improve Self Image and Reduce Peritonitis*



*MY-733 and MY-744 are the economy versions of MY-730 and MY-740 respectively. They both use the solution bag from the previous exchange to serve as the drain bag. Sterile connectors are supplied with the package to cap off the empty solution bag.*

*MY-730 is designed for use with solution bags which have a spike type input port, while MY-740 is for use with solution bags which have a luer type input port.*

### SAFETY

- Flush-Before-Fill reduces peritonitis
- Pre-filled disinfectant cap (QC Cap) ensures dependable aseptic connection
- Recessed and shielded connectors prevent touch contamination
- Colour coded clamps & connectors minimize operational errors
- Closed system eliminates open drain containers and minimizes contamination
- Non-wearable set increases patient compliance

### CONVENIENCE

- No bulky bags and no long tubing sets to wear
- Pre-filled disinfectant cap (QC Cap) provides quick patient connect and disconnect
- Improves patient's lifestyle and self-image
- Does not require a prep-kit
- No external disinfectant required
- No droppers or dispensers required
- No extra caps required

