

# ***EASY CARE TUBING SETS*** ***and accessories***

**Easy Care Tubing Set**



## **EASY CARE TUBING SET**

- Designed for use with Easy Care Peritoneal dialysis machine
- Drop-in auto-lock cassette simplifies tubing set loading and operation
- Integrated drainage bag provides closed loop system to minimize contamination
- QC cap system ensures sterility during and after treatment

**QC CAP QC-200**



- QC Cap Provides quick patient connection and disconnection
- Pre-filled disinfectant cap ensures dependable aseptic connection
- Provides continuous disinfection of patient connector between dialysis treatments
- Recessed connector prevents touch contamination
- Reduces peritonitis

## **TRANSFER SET QT-012C**



Medionics' transfer set connects to all catheter adapters, and is an essential component of the Medionics' QC System.



MEDIONICS INTERNATIONAL INC.  
114 Anderson Avenue, Markham, Ontario L6E 1A5 Canada  
Tel: (905) 472-6544 [info@medionics.com](mailto:info@medionics.com) [www.medionics.com](http://www.medionics.com)

